

# LZI Oral Fluid Ecstasy (MDMA) Calibrators



**FOR RESEARCH & DEVELOPMENT USE ONLY**

## ***Lin-Zhi International, Inc.***

<b>REF</b>	<b>Description</b>	<b>Quantity</b>
S0001	LZI Oral Fluid Negative Calibrator	1 x 5 mL
S0167	LZI Oral Fluid Ecstasy (MDMA) 15 ng/mL Low Calibrator / Level 1 Control	1 x 5 mL
S0163	LZI Oral Fluid Ecstasy (MDMA) 25 ng/mL Cutoff Calibrator	1 x 5 mL
S0168	LZI Oral Fluid Ecstasy (MDMA) 35 ng/mL Intermediate Calibrator / Level 2 Control	1 x 5 mL
S0165	LZI Oral Fluid Ecstasy (MDMA) 50 ng/mL High Calibrator	1 x 5 mL

### **Intended Use**

The Lin-Zhi International, Inc. (LZI) Oral Fluid Ecstasy (MDMA) Calibrators are for use as calibrators in the qualitative and semi-quantitative calibration of the LZI Oral Fluid Ecstasy (MDMA) Enzyme Immunoassay (Ref# S0160/S0161) on a number of automated clinical chemistry analyzers.

### **Description of the Calibrators:**

The LZI Oral Fluid Ecstasy (MDMA) Calibrators are in a negative synthetic oral fluid matrix, and ready to use. The LZI Oral Fluid Negative Calibrator is a drug-free synthetic oral fluid matrix containing buffers, stabilizers, and less than 0.1 % of sodium azide. The calibrators are prepared by spiking known concentrations of 3, 4-methylenedioxymethamphetamine into the LZI Oral Fluid Negative Calibrator. Calibrators are made from NIST traceable standards.

\*Actual concentrations of these calibrators are within  $\pm 10\%$  of the target value as determined by GC/MS. Values are provided only as guidelines and laboratories should determine the ranges based on their own test system and tolerance (1).

### **Precautions and Warnings**

- *The LZI Oral Fluid Ecstasy Calibrators are Non-FDA approved and are for Research & Development Use Only. The LZI Oral Fluid Ecstasy Calibrators should not be repackaged for in vitro diagnostic use.*
- *Harmful if swallowed.*
- *The calibrators contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such liquids always flush with a large volume of water to prevent azide build-up.*
- *The calibrators are prepared from a non-sterile synthetic oral fluid matrix. Always apply good laboratory practices to avoid any skin contact or ingestion.*
- *Do not use the calibrators beyond their expiration dates.*

### **Preparation and Storage**

The calibrators are ready-to-use. No reconstitution is required. Label the cap before removal to identify it with the original bottle. The calibrators should be stored refrigerated at 2-8 °C when not in use.

### **Stability**

When stored refrigerated at 2-8 °C, the calibrators are stable either opened and recapped, or unopened, until the expiration date printed on the vial label. Store calibrators tightly capped when not in use. Calibrator solution dispensed in the sample cups and left on board the clinical analyzer should be discarded after use.

### **Procedure and Results**

For qualitative calibration, use the 25 ng/mL as the cutoff calibrator. For semi-quantitative calibration, use all five calibrators. Recalibration should be performed after reagent bottle change, a change in calibrators or reagent lot, and after instrument maintenance is performed. For interpretation of results, refer to the appropriate LZI Oral Fluid Ecstasy (MDMA) Enzyme Immunoassay (Ref# S0160/S0161) package insert (2).

### **Limitations**

The LZI Oral Fluid Ecstasy (MDMA) Calibrators are for use with the LZI Oral Fluid Ecstasy (MDMA) Enzyme Immunoassay (Ref# S0160/S0161) to detect 3, 4-methylenedioxymethamphetamine in human oral fluid only.

### **Bibliography**

1. Guidance for Industry, Abbreviated 510(k) Submissions for In Vitro Diagnostic Calibrators. U.S. Department of Health and Human Services. FDA, Document issued on February 22, 1999.
2. LZI Oral Fluid Ecstasy (MDMA) Enzyme Immunoassay (Ref# S0160/S0161) package insert.

**Notice: Adulteration of reagents, use of instruments without appropriate capabilities, or other failure to follow instructions as set forth in this labeling can affect performance characteristics and stated or implied label claims.**



### **Manufacturer:**

**Lin-Zhi International, Inc.**  
2945 Oakmead Village Court  
Santa Clara, CA 95051  
USA  
Tel: (408) 970-8811  
Fax: (408) 970-9030  
www.lin-zhi.com