

# LZI Ethyl Glucuronide Calibrators



| For Forensic Use Only

## Lin-Zhi International, Inc.

REF	Description	Quantity
0001	Negative Calibrator	1 x 5 mL
0322b	Ethyl Glucuronide 250 ng/mL Low Calibrator	1 x 5 mL
0323b	Ethyl Glucuronide 1000 ng/mL Cutoff Calibrator	1 x 5 mL
0324b	Ethyl Glucuronide 2000 ng/mL Intermediate Calibrator #1	1 x 5 mL
0325b	Ethyl Glucuronide 3000 ng/mL Intermediate Calibrator #2	1 x 5 mL
0326b	Ethyl Glucuronide 4000 ng/mL High Calibrator	1 x 5 mL

### Intended Use

The Lin-Zhi International, Inc. (LZI) Ethyl Glucuronide Calibrators are for use as calibrators in the qualitative and semi-quantitative calibration of the LZI Ethyl Glucuronide Enzyme Immunoassay (Ref# 0320b/0321b) on a number of automated clinical chemistry analyzers (1). These calibrators are for Forensic Use Only and should not be repackaged for *in vitro* diagnostic use.

### Description of the Calibrators:

The LZI Ethyl Glucuronide Calibrators are human urine-based liquids, and ready-to-use. The Negative Calibrator is a processed drug-free human urine matrix containing buffers, stabilizers, and less than 0.1 % of sodium azide. The calibrators are prepared by spiking known concentrations of ethyl glucuronide into the Negative Calibrator.

### Precautions and Warning

- The LZI Ethyl Glucuronide Calibrators should not be repackaged for *in vitro* diagnostic use.
- Harmful if swallowed.
- The calibrators contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azide. When disposing such liquids always flush with a large volume of water to prevent azide build-up (2).
- The calibrators are prepared from non-sterile human urine. They are not tested by licensed reagents for the presence of antibodies to human immunodeficiency viruses, the hepatitis antigens, and/or anti-hepatitis antibodies. They should be handled as potentially infectious. Always apply good laboratory practice to avoid any skin contact or ingestion.
- Do not use the calibrators beyond their expiration dates.

### Preparation and Storage

The calibrators are provided ready-to-use. No reconstitution is required. Label the cap before removal to identify it with the original bottle. The calibrators should be stored refrigerated at 2-8°C when not in use.

### Stability

When stored refrigerated at 2-8°C, the calibrators are stable either opened-recapped or unopened until the expiration date printed on the vial label. Store calibrators tightly capped when not in use. Calibrator solution dispensed in the sample cups and left on board of the clinical analyzer should be discarded after use.

### Procedure and Results

For qualitative calibration, use the 1000 ng/mL as your cutoff calibrator. For semi-quantitative calibration, use all six calibrators. Recalibration should be performed after reagent bottle change or there is a change in calibrators or reagent lot, and after instrument maintenance is performed. For interpretation of results, refer to the appropriate LZI Ethyl Glucuronide Enzyme Immunoassay (Ref# 0320b/0321b) package insert (1).

### Limitations

The LZI Ethyl Glucuronide Calibrators are for use with the LZI Ethyl Glucuronide Enzyme Immunoassay (Ref# 0320b/0321b) to detect ethyl glucuronide in human urine only.

### Bibliography

1. LZI Ethyl Glucuronide Enzyme Immunoassay (Ref# 0320b/0321b) package insert.
2. Sodium Azide. National Institute for Occupational Safety (NIOSH). Pocket Guide to Chemical Hazards. Third Printing, September 2007. Available online at: <https://www.cdc.gov/niosh/npg/default.html>

**Notice: Adulteration of reagents, use of instruments without appropriate capabilities, or other failure to follow instructions as set forth in this labeling can affect performance characteristics, and stated or implied label claims.**

Additions, deletions, or changes are indicated by a change bar in the margin.



### Manufacturer:

Lin-Zhi International, Inc.  
2945 Oakmead Village Court  
Santa Clara, CA 95051  
USA  
Tel: (408) 970-8811  
Fax: (408) 970-9030  
[www.lin-zhi.com](http://www.lin-zhi.com)