

LZI Methadone Enzyme Immunoassay

IVD For In Vitro Diagnostic Use Only



REF 0110 (100/37.5 mL R₁/R₂ Kit)
0111 (1000/375 mL R₁/R₂ Kit)



Lin-Zhi International, Inc.

Intended Use

The Lin-Zhi International, Inc. (LZI) Methadone Enzyme Immunoassay is intended for the qualitative and semi-quantitative determination of methadone in human urine at a cutoff value of 300 ng/mL. The assay is designed for prescription use with a number of automated clinical chemistry analyzers.

The assay provides only a preliminary analytical result. A more specific alternative analytical chemistry method must be used in order to obtain a confirmed analytical result. Gas or Liquid Chromatography/Mass Spectrometry (GC/MS or LC/MS) are the preferred confirmatory methods (1, 2). Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.

Summary and Explanation of Test

Methadone is a synthetic diphenylheptanonylamine opioid that has similar analgesic activity and potency to morphine when administered parenterally. However, unlike morphine, it reliably retains its effectiveness when given orally, and tolerance and physical dependency develop slowly (3, 4). Although methadone is prescribed to relieve chronic pain, its primary application is the detoxification and/or treatment of narcotic or heroin addiction (3-6). The abuse potential of methadone is comparable to that of morphine due to its similar pharmacological activity (3, 5, 7).

Methadone is available in tablets and as a solution for parenteral injection. It is readily absorbed from the gastrointestinal tract when ingested, and metabolized extensively in the liver. Initial N-demethylation results in normethadone, which rapidly undergoes cyclization followed by dehydration to form the 2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine, commonly known as EDDP. Further N-demethylation yields a secondary metabolite, the 2-ethyl-5-methyl-3,3-diphenyl-1-pyrroline (EMDP) (8). The metabolites are secreted in urine or bile along with unchanged drug.

Assay Principle

The LZI Methadone assay is a homogeneous enzyme immunoassay ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent (9). Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity.

In the absence of drug in the sample, methadone-labeled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when drug is present in the sample, antibody binds to free drug; the unbound methadone-labeled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at a 340 nm primary wavelength.

Reagents Provided

Antibody/Substrate Reagent (R₁): Contains mouse monoclonal anti-methadone antibody, glucose-6-phosphate (G6P), nicotinamide adenine dinucleotide (NAD), stabilizers, and sodium azide (0.09 %) as a preservative.

Enzyme-drug Conjugate Reagent (R₂): Contains glucose-6-phosphate dehydrogenase (G6PDH) labeled with methadone in buffer with sodium azide (0.09 %) as a preservative.

Calibrators and Controls*

*Calibrators and controls are sold separately and contain negative human urine with sodium azide as a preservative.

| METHADONE Calibrators | REF |
|---|------|
| Negative Calibrator | 0001 |
| Low Calibrator: Contains 150 ng/mL methadone | 0112 |
| Cutoff Calibrator: Contains 300 ng/mL methadone | 0113 |
| Intermediate Calibrator: Contains 600 ng/mL methadone | 0114 |
| High Calibrator: Contains 1000 ng/mL methadone | 0115 |
| METHADONE Controls | REF |
| Level 1 Control: Contains 225 ng/mL methadone | 0117 |
| Level 2 Control: Contains 375 ng/mL methadone | 0118 |

Precautions and Warning

- This test is for in vitro diagnostic use only. Harmful if swallowed.
- Reagent contains sodium azide as a preservative, which may form explosive compounds in metal drain lines. When disposing such reagents or wastes always flush with a large volume of water to prevent azide build-up. See National Institute for Occupational Safety and Health Bulletin: Explosive Azide Hazards (10).
- Do not use the reagents beyond their expiration dates.
- For USA: Caution: Federal law restricts this device to sale by or on the order of a physician.

Reagent Preparation and Storage

The reagents are ready-to-use. No reagent preparation is required. All assay components should be refrigerated at 2-8°C when not in use.

Specimen Collection and Handling

Urine samples may be collected in plastic or glass containers. Some plastics may absorb drugs. Use of plastics such as polyethylene is recommended (11). Use fresh urine specimens for the test. If a sample cannot be analyzed immediately, it may be refrigerated at 2-8°C for up to seven days (12). For longer storage, keep sample frozen at -20°C and then thaw before use. Studies have shown methadone analytes in urine are stable at -20°C up to 384 days (13). Samples should be at room temperature (18-25°C) for testing. Samples with high turbidity should be centrifuged before analysis. Adulteration may cause erroneous results. If sample adulteration is suspected, obtain a new sample and forward both samples to the laboratory for testing. Handle all urine specimens as if they are potentially infectious.

Instrument

Clinical chemistry analyzers capable of maintaining a constant temperature, pipetting samples, mixing reagents, measuring enzyme rates at a 340 nm primary wavelength and timing the reaction accurately can be used to perform this homogeneous immunoassay. Performance characteristics presented in this package insert have been validated on the Synchron® CX4CE and the Synermed IR500 clinical analyzers.

Assay Procedure

Refer to the specific parameters used for each analyzer before performing the assay. For qualitative analysis, use the 300 ng/mL as the cutoff calibrator. For semi-quantitative analysis, use all five calibrators. Recalibration should be performed after reagent bottle change or if there is a change in calibrators or reagent lot. Two levels of controls are also available for monitoring the cutoff level: 225 ng/mL and 375 ng/mL.

Calibration and Quality Control

Good laboratory practices recommend the use of at least two levels of control specimens (one positive and one negative control near the cutoff) to ensure proper assay performance. Controls should be run with each new calibration and after specific maintenance or troubleshooting procedures as detailed in the instrument system manual. Each laboratory should establish its own control frequency. If any trends or sudden change in control value are observed, review all operating parameters, or contact LZI technical support for further assistance. Laboratories should comply with all federal, state, and local laws, as well as all guidelines and regulations.

Results

Note: A preliminary positive test result does not necessarily mean a person took illegal drugs and a negative test result does not necessarily mean a person did not take illegal drugs. There are a number of factors that influence the reliability of drug tests.

Qualitative: The cutoff calibrator which contains 300 ng/mL of methadone is used as a reference for distinguishing a preliminary positive from negative samples. A sample with a change in absorbance ($\Delta A/\text{min}$) equal to or greater than that obtained with the cutoff calibrator is considered a preliminary positive. A sample with a change in absorbance ($\Delta A/\text{min}$) lower than that obtained with the cutoff calibrator is considered negative.

Semi-Quantitative: The semi-quantitative mode is for purposes of (1) enabling laboratories to determine an appropriate dilution of the specimen for verification by a confirmatory method such as GC/MS, LC/MS or (2) permitting laboratories to establish quality control procedures.

When an approximation of concentration is required, a calibration curve can be established with five calibrators. The concentration of methadone in the sample may then be estimated from the calibration curve.

Limitations

1. A preliminary positive result from the assay indicates only the presence of methadone. The test is not intended for quantifying this single analyte in samples.
2. A preliminary positive result does not necessarily indicate drug abuse.
3. A negative result does not necessarily mean a person did not take illegal drugs.
4. Care should be taken when reporting results as numerous factors (e.g., fluid intake, endogenous or exogenous interferents) may influence the urine test result.
5. Preliminary positive results should be confirmed by other affirmative, analytical chemistry methods (e.g., chromatography), preferably GC/MS or LC/MS.
6. The test is designed for use with human urine only.
7. The test is not for therapeutic drug monitoring.

Typical Performance Characteristics

The results shown below were obtained with the Synchron CX4CE clinical chemistry analyzer.

Precision:

Qualitative analysis: The three calibrators and two levels of controls were evaluated. Typical results (Δ mA/min) are as follows:

| Concentration | Within Run (N=21) | | | Run-to-Run* (N=12) | | |
|---------------|-------------------|-----|-------|--------------------|-----|-------|
| | Mean | SD | % CV | Mean | SD | % CV |
| 0 ng/mL | 209.4 | 1.0 | 0.5 % | 209.5 | 1.1 | 0.5 % |
| 225 ng/mL | 272.8 | 1.1 | 0.4 % | 271.4 | 2.0 | 0.7 % |
| 300 ng/mL | 293.4 | 0.7 | 0.3 % | 292.3 | 1.9 | 0.7 % |
| 375 ng/mL | 308.2 | 0.9 | 0.3 % | 307.2 | 2.6 | 0.8 % |
| 1000 ng/mL | 344.9 | 1.1 | 0.3 % | 344.1 | 2.0 | 0.6 % |

*Run-to-Run testing completed over 3 weeks

Semi-quantitative analysis: The concentrations of the cutoff level and the two levels of controls were determined with reference curves from five calibrators. Typical results (ng/mL) are as follows:

| Concentration | Within Run (N=21) | | | Run-to-Run* (N=12) | | |
|---------------|-------------------|-----|-------|--------------------|------|-------|
| | Mean | SD | % CV | Mean | SD | % CV |
| 225 ng/mL | 229.2 | 2.5 | 1.1 % | 228.0 | 7.8 | 3.4 % |
| 300 ng/mL | 308.5 | 3.5 | 1.1 % | 305.1 | 10.2 | 3.4 % |
| 375 ng/mL | 391.4 | 5.4 | 1.4 % | 379.6 | 13.2 | 3.5 % |

*Run-to-Run testing completed over 3 weeks

Sensitivity: Sensitivity, defined as the lowest concentration that can be differentiated from negative urine with 95 % confidence, was tested to be 15 ng/mL.

Accuracy: One hundred and sixty (160) clinical urine specimens were tested with current EIA, 49 samples were positive, and 111 samples were negative. All positive samples were confirmed with GC/MS results.

| Cutoff Value (300 ng/mL) | GC/MS | LZI MTD EIA | % Agreement with Predicate |
|--------------------------|-------|-------------|----------------------------|
| # Positive Samples | 49 | 49 | 100 % |
| # Negative Samples | 111 | 111 | 100 % |
| Total # of Samples | 160 | 160 | N/A |

In addition to the above study, 19 diluted clinical samples with a methadone concentration ranging from 205 ng/mL to 434 ng/mL (determined by GC/MS) were evaluated with the LZI Methadone EIA. The eight samples with methadone GC/MS values greater than the cutoff (ranging from 305 to 434 ng/mL) tested positive by EIA. Among the 11 samples with methadone GC/MS below the cutoff (ranging from 205 to 288 ng/mL), three samples were found negative while eight samples tested positive by the EIA. The GC/MS concentrations of methadone in these eight samples range from 214 ng/mL to 288 ng/mL. Each of these samples also contained a substantial amount of the primary methadone metabolite EDDP (ranging from 247 ng/mL to 825 ng/mL).

Analytical Recovery: Analytical recovery was evaluated by spiking known concentrations of methadone to negative urine samples.

In qualitative analysis, the assay correctly identified spiked samples containing more than 300 ng/mL of methadone (n=25, spiked levels equal to or higher than Level 2 Control, 375 ng/mL) as positive, and those containing less than 300 ng/mL of methadone (n=25, spiked levels equal to or less than Level 1 Control, 225 ng/mL) as negative.

For semi-quantitative analysis, the average recovery for samples spiked with 30 to 900 ng/mL of methadone (five samples at each level) is summarized in the following table:

Analytical Recovery, continued:

| Expected Value (ng/mL) | Observed Value (ng/mL) | % Recovery |
|------------------------|------------------------|------------|
| 30 | 29.5 | 98.5 % |
| 60 | 64.2 | 106.9 % |
| 120 | 122.9 | 102.5 % |
| 180 | 176.3 | 97.9 % |
| 225 | 217.3 | 96.6 % |
| 375 | 351.5 | 93.7 % |
| 400 | 460.5 | 92.1 % |
| 600 | 531.3 | 88.5 % |
| 750 | 638.6 | 85.1 % |
| 900 | 700.8 | 77.9 % |

Specificity: Various potentially interfering substances were tested for cross-reactivity with the assay. Test compounds were spiked into the drug-free urine calibrator matrix to various concentrations and evaluated against the cutoff calibrator.

The table below lists the concentration of each test compound that gave a response approximately equivalent to that of the cutoff calibrator (as positive) or the maximal concentration of the compound tested that gave a response below the response of the cutoff calibrator (as negative).

Structurally Related Methadone Compounds:

| Compound | Target [μ g/mL] | % Cross-Reactivity |
|----------------------------|----------------------|--------------------|
| Methadone | 0.30 | Positive |
| LAAM-HCl | 10 | Positive |
| (-) α -Methadol-HCl | 8 | Positive |
| nor-LAAM-HCl | 10 | Negative |
| EDDP-HI | 100 | Negative |
| EMDP-HCl | 100 | Negative |

Structurally Unrelated Pharmacological Compounds:

| Compound | Target [μ g/mL] | % Cross-Reactivity |
|----------------------|----------------------|--------------------|
| Acetaminophen | 1000 | Negative |
| Acetylsalicylic Acid | 1000 | Negative |
| Amitriptyline | 1000 | Negative |
| Amobarbital | 1000 | Negative |
| Amphetamine | 1000 | Negative |
| Benzoyllecgonine | 1000 | Negative |
| Bupropion | 1000 | Negative |
| Caffeine | 1000 | Negative |
| Chlorpheniramine | 1000 | Negative |
| Chlorpromazine | 1000 | Negative |
| Cocaine | 1000 | Negative |
| Codeine | 1000 | Negative |
| Dextromethorphan | 1000 | Negative |
| Egonine | 1000 | Negative |
| Ephedrine | 1000 | Negative |
| Ibuprofen | 2000 | Negative |
| Imipramine | 1000 | Negative |
| Lidocaine | 1000 | Negative |
| Meperidine | 1000 | Negative |
| Methamphetamine | 1000 | Negative |
| Methaqualone | 1000 | Negative |
| Morphine | 1000 | Negative |
| Nortriptyline | 1000 | Negative |
| Oxazepam | 1000 | Negative |
| Phencyclidine | 1000 | Negative |
| Phenobarbital | 1000 | Negative |
| Promethazine | 1000 | Negative |
| Propoxyphene | 1000 | Negative |
| Ranitidine | 1000 | Negative |
| Secobarbital | 1000 | Negative |
| Valproic Acid | 1000 | Negative |

It is possible that other substances and/or factors not listed above may interfere with the test and cause false positive results.

Bibliography:

1. Urine Testing for Drug of Abuse, National Institute on Drug Abuse (NIDA) Research Monograph 73, (1986).
2. Mandatory Guidelines for Federal Workplace Drug Testing Program, National Institute on Drug Abuse, Federal Register, 53(69): 11970 (1988).
3. Goodman, L.S., and Gilman, A., The Pharmacological Basis of Therapeutics, 4th edition, 380, The MacMillan Co., (1970).
4. Roper, J.D., Garside, D., and Goldberger, B.A., "Opiates", in Contemporary Practice in Clinical Toxicology, 2nd edition, Leslie M. Shaw, editor-in-chief, AACCC, (2000).
5. Lacy, C., Armstrong, L.L., Lipsy, R.J., and Lance, L.L., Drug Information Handbook. Hudson, OH: Lexi-Comp, (1993).
6. Katzung, B.G., Basic and Clinical Pharmacology, 6th edition, Appleton & Lange, Norwalk, CT. (1995).

Bibliography, continued:

7. Physician's Desk Reference, 54th edition, Medical Economics Company, Montvale, NJ, 2711-2713 (2000).
8. Baselt, R.C., and Carvey, R.H., Disposition of Toxic Drugs and Chemicals in Man, 4th edition, Chemical Toxicology Institute, Foster City, CA. (1995).
9. Rubenstein, K.E., Schneider, R.S., and Ullman, E.F., Homogeneous Enzyme Immunoassay: A New Immunochemical Technique, *Biochem Biophys Res Commun*, **47**:846 (1972).
10. Sodium Azide. National Institute for Occupational Safety (NIOSH). Pocket Guide to Chemical Hazards. Third Printing, September 2007. Available online at: <https://www.cdc.gov/niosh/npg/default.html>
11. Yahya, A.M., McElmay, J.C., and D'Arcy, P.F., Drug absorption to glass and plastics, *Drug Metabol Drug Interact*, **6**(1):1-45 (1988).
12. Ciuti, R., Quercioli, M., and Borsotti, M., Drug of abuse stability in native urine specimens vs. stabilized urine samples, *Biochimica Clinica*, **38**(2) 103-109 (2014). Moody, D.E., Monti, K.M., and Spanbauer, A.C., Long-Term Stability of Abused Drugs and Antiabuse Chemotherapeutical Agents Stored at -20°C, *J. Anal. Toxicol.*, **23**:535-540 (1999).
13. Moody, D.E., Monti, K.M., Spanbauer, A.C., and Hsu, J.P., Long-Term Stability of Abused Drugs and Antiabuse Chemotherapeutical Agents Stored at -20°C, *J Anal Toxicol*. **23**:535-540 (1999).

Additions, deletions, or changes are indicated by a change bar in the margin.
For technical assistance please call: (408) 970-8811

Manufacturer:

Lin-Zhi International, Inc.
2945 Oakmead Village Court
Santa Clara, CA 95051
USA
Tel: (408) 970-8811
Fax: (408) 970-9030
www.lin-zhi.com



**Authorized European
Rep. within the EU:**

CEpartner4U
Esdoornlaan 13
3951DB Maarn
The Netherlands
www.cepartner4u.eu



© June 2019 Rev. 13

Printed in USA