



PRODUCT CATALOG: Drugs of Abuse Oral Fluid by Alphabetical Order

| <u>REF #</u> | <u>Product Name</u> | <u>REF #</u> | <u>Product Name</u> |
|--------------|--|--------------|--|
| S0000b | LZI Oral Fluid Collector, 500 tubes/case | S0220 | Ethyl Alcohol Assay Kit, small test kit |
| S0002 | Oral Fluid Recovery Buffer, 500 or 1000 mL | S0221 | Ethyl Alcohol Assay Kit, large test kit |
| S0005 | Oral Fluid 100 Sarstedt Collector Kit | S0001 | Oral Fluid Negative Calibrator |
| S0006 | Oral Fluid 1000 Sarstedt Collector Kit | S0223 | Ethyl Alcohol Assay Calibrator (50 mg/dL) |
| S0290 | 6-Acetylmorphine EIA Kit, small test kit | S0226 | Ethyl Alcohol Assay Calibrator/Control (25 mg/dL) |
| S0291 | 6-Acetylmorphine EIA Kit, large test kit | S0227 | Ethyl Alcohol Assay Calibrator/Control (75 mg/dL) |
| S0298 | 6-Acetylmorphine Negative Calibrator (0 ng/mL) | S0228 | Ethyl Alcohol Assay Calibrator (100 mg/dL) |
| S0292 | 6-Acetylmorphine Low Calibrator/Level 1 Control (2 ng/mL) | S0110c | Methadone EIA Kit, small test kit |
| S0293 | 6-Acetylmorphine Cutoff Calibrator (4 ng/mL) | S0111c | Methadone EIA Kit, large test kit |
| S0294 | 6-Acetylmorphine Intermediate Calibrator (10 ng/mL) | S0008 | Oral Fluid LIS Negative Calibrator |
| S0295 | 6-Acetylmorphine High Calibrator (20 ng/mL) | S0113c | Methadone Cutoff Calibrator (20 ng/mL) |
| S0297 | 6-Acetylmorphine Level 2 Control (6 ng/mL) | S0115c | Methadone High Calibrator (50 ng/mL) |
| S0040b | ^{510(k)} _{Clear} Amphetamine EIA Kit, small test kit | S0117c | Methadone Level 1 Control/Low Calibrator (10 ng/mL) |
| S0041b | Amphetamine EIA Kit, large test kit | S0118c | Methadone Level 2 Control/Intermediate Calibrator (30 ng/mL) |
| S0001 | Oral Fluid Negative Calibrator | S0190 | Methadone Metabolite (EDDP) EIA Kit, small test kit |
| S0042b | Amphetamine Low Calibrator (20 ng/mL) | S0191 | Methadone Metabolite (EDDP) EIA Kit, large test kit |
| S0043b | Amphetamine Cutoff Calibrator (50 ng/mL) | S0001 | Oral Fluid Negative Calibrator |
| S0044b | Amphetamine Intermediate Calibrator (100 ng/mL) | S0193 | Methadone Metabolite Cutoff Calibrator (20 ng/mL) |
| S0045b | Amphetamine High Calibrator (140 ng/mL) | S0195 | Methadone Metabolite High Calibrator (50 ng/mL) |
| S0046b | Amphetamine Level 1 Control (37.5 ng/mL) | S0197 | Methadone Metabolite Level 1 Control/Low Calibrator (10 ng/mL) |
| S0047b | Amphetamine Level 2 Control (62.5 ng/mL) | S0198 | Methadone Metabolite Level 2 Control/Intermediate Cal. (30 ng/mL) |
| S0140b | Barbiturate EIA Kit, small test kit | S0050b | ^{510(k)} _{Clear} Methamphetamine EIA Kit, small test kit |
| S0141b | Barbiturate EIA Kit, large test kit | S0051b | Methamphetamine EIA Kit, large test kit |
| S0001 | Oral Fluid Negative Calibrator | S0001 | Oral Fluid Negative Calibrator |
| S0143b | Barbiturate Cutoff Calibrator (Secobarbital 20 ng/mL) | S0052b | Methamphetamine Low Calibrator (20 ng/mL) |
| S0145b | Barbiturate High Calibrator (Secobarbital 50 ng/mL) | S0053b | Methamphetamine Cutoff Calibrator (50 ng/mL) |
| S0147b | Barbiturate Level 1 Control/Low Cal. (Secobarbital 10 ng/mL) | S0054b | Methamphetamine Intermediate Calibrator (100 ng/mL) |
| S0148b | Barbiturate Level 2 Control/Inter. Cal. (Secobarbital 30 ng/mL) | S0055b | Methamphetamine High Calibrator (140 ng/mL) |
| S0130c | Benzodiazepine EIA Kit, small test kit | S0056b | Methamphetamine Level 1 Control (37.5 ng/mL) |
| S0131c | Benzodiazepine EIA Kit, large test kit | S0057b | Methamphetamine Level 2 Control (62.5 ng/mL) |
| S0008 | Oral Fluid LIS Negative Calibrator | S0020c | Opiate EIA Kit, small test kit |
| S0133c | Benzodiazepine Cutoff Calibrator (Oxazepam 20 ng/mL) | S0021c | Opiate EIA Kit, large test kit |
| S0135c | Benzodiazepine High Calibrator (Oxazepam 50 ng/mL) | S0008 | Oral Fluid LIS Negative Calibrator |
| S0137c | Benzodiazepine Level 1 Control/Low Cal. (Oxazepam 10 ng/mL) | S0022c | Opiate Low Calibrator (Morphine 20 ng/mL) |
| S0138c | Benzodiazepine Level 2 Control/Inter. Cal. (Oxazepam 30 ng/mL) | S0023c | Opiate Cutoff Calibrator (Morphine 40 ng/mL) |
| S0030d | Cocaine Metabolite EIA Kit, small test kit | S0024c | Opiate Intermediate Calibrator (Morphine 60 ng/mL) |
| S0031d | Cocaine Metabolite EIA Kit, large test kit | S0025c | Opiate High Calibrator (Morphine 80 ng/mL) |
| S0007d | Oral Fluid LC Negative Calibrator | S0027c | Opiate Level 1 Control (Morphine 30 ng/mL) |
| S0033d | Cocaine Metabolite Low Calibrator (BE 10 ng/mL) | S0028c | Opiate Level 2 Control (Morphine 50 ng/mL) |
| S0034d | Cocaine Metabolite Cutoff Calibrator (BE 20 ng/mL) | S0240b | Oxycodone EIA Kit, small test kit |
| S0035d | Cocaine Metabolite Intermediate Calibrator (BE 35 ng/mL) | S0241b | Oxycodone EIA Kit, large test kit |
| S0036d | Cocaine Metabolite High Calibrator (BE 50 ng/mL) | S0001 | Oral Fluid Negative Calibrator |
| S0037d | Cocaine Metabolite Level 1 Control (BE 15 ng/mL) | S0242b | Oxycodone Low Calibrator (20 ng/mL) |
| S0038d | Cocaine Metabolite Level 2 Control (BE 25 ng/mL) | S0243b | Oxycodone Cutoff Calibrator (40 ng/mL) |
| S0230 | Cotinine EIA Kit, small test kit | S0244b | Oxycodone Intermediate Calibrator (60 ng/mL) |
| S0231 | Cotinine EIA Kit, large test kit | S0245b | Oxycodone High Calibrator (80 ng/mL) |
| S0001 | Oral Fluid Negative Calibrator | S0247b | Oxycodone Level 1 Control (30 ng/mL) |
| S0233 | Cotinine Cutoff Calibrator (25 ng/mL) | S0248b | Oxycodone Level 2 Control (50 ng/mL) |
| S0235 | Cotinine High Calibrator (75 ng/mL) | S0010c | Phencyclidine EIA Kit, small test kit |
| S0237 | Cotinine Level 1 Control/Low Calibrator (10 ng/mL) | S0011c | Phencyclidine EIA Kit, large test kit |
| S0238 | Cotinine Level 2 Control/Intermediate Calibrator (50 ng/mL) | S0008 | Oral Fluid LIS Negative Calibrator |
| S0160 | Ecstasy (MDMA) EIA Kit, small test kit | S0012c | Phencyclidine Low Calibrator (5 ng/mL) |
| S0161 | Ecstasy (MDMA) EIA Kit, large test kit | S0013c | Phencyclidine Cutoff Calibrator (10 ng/mL) |
| S0001 | Oral Fluid Negative Calibrator | S0014c | Phencyclidine Intermediate Calibrator (15 ng/mL) |
| S0163 | Ecstasy Cutoff Calibrator (25 ng/mL) | S0015c | Phencyclidine High Calibrator (25 ng/mL) |
| S0165 | Ecstasy High Calibrator (50 ng/mL) | S0017c | Phencyclidine Level 1 Control (7.5 ng/mL) |
| S0167 | Ecstasy Level 1 Control/Low Calibrator (15 ng/mL) | S0018c | Phencyclidine Level 2 Control (12.5 ng/mL) |
| S0168 | Ecstasy Level 2 Control/Intermediate Calibrator (35 ng/mL) | | |

^{510(k)}_{Clear} = Received Substantially Equivalent 510(k) decision from FDA

* All Oral Fluid kit formulations are for Forensic Use Only (Non-FDA Approved) unless denoted otherwise with an ^{510(k)}_{Clear} icon.

* The Negative Calibrator to be used with each specific assay has been listed with each assay set.

* All single analyte calibrators and controls are 5 mL per bottle.


* A small test kit is 75/37.5 mL R₁/R₂


* A large test kit is 750/375 mL R₁/R₂



PRODUCT CATALOG: Drugs of Abuse Oral Fluid by Alphabetical Order

| <u>REF #</u> | <u>Product Name</u> |
|--------------|--|
| S0120c | Propoxyphene EIA Kit, small test kit |
| S0121c | Propoxyphene EIA Kit, large test kit |
| S0008 | Oral Fluid LIS Negative Calibrator |
| S0123c | Propoxyphene Cutoff Calibrator (20 ng/mL) |
| S0125c | Propoxyphene High Calibrator (50 ng/mL) |
| S0127c | Propoxyphene Level 1 Control/Low Calibrator (10 ng/mL) |
| S0128c | Propoxyphene Level 2 Control/Intermediate Calibrator (30 ng/mL) |
| S0070 | THC (Cannabinoid) EIA Kit, small test kit |
| S0071 | THC (Cannabinoid) EIA Kit, large test kit |
| S0072 | Oral Fluid Cannabinoid Negative Calibrator |
| S0073 | THC Cutoff Calibrator (Δ^9 -THC 4 ng/mL) |
| S0075 | THC High Calibrator (Δ^9 -THC 20 ng/mL) |
| S0077 | THC Level 1 Control/Low Calibrator (Δ^9 -THC 2 ng/mL) |
| S0078 | THC Level 2 Control/Intermediate Calibrator (Δ^9 -THC 8 ng/mL) |

 = Received Substantially Equivalent 510(k) decision from FDA

* All Oral Fluid kit formulations are for Forensic Use Only (Non-FDA Approved) unless denoted otherwise with an  icon.

* The Negative Calibrator to be used with each specific assay has been listed with each assay set.

* All single analyte calibrators and controls are 5 mL per bottle.

* A small test kit is 75/37.5 mL R₁/R₂

* A large test kit is 750/375 mL R₁/R₂