

Multi-Analyte Calibrators and Controls Set A

<i>Works with 5 EIAs</i>	#0812 Low Calibrator	#0817 Level I Control	#0813 Cutoff Calibrator	#0818 Level II Control	#0814 Intermediate Calibrator	#0815 High Calibrator
Benzoyllecgonine (Cocaine 150 EIA)		112.5 ng/mL	150 ng/mL	187.5 ng/mL	300 ng/mL	1000 ng/mL
d-Methamphetamine (Amphetamines 500 EIA & Methamphetamine EIA)	250 ng/mL	375 ng/mL	500 ng/mL	625 ng/mL	750 ng/mL	1000 ng/mL
Methadone (Methadone EIA)	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
6-Acetylmorphine (6-Acetylmorphine EIA)	5 ng/mL	7.5 ng/mL	10 ng/mL	12.5 ng/mL	20 ng/mL	40 ng/mL

Multi-Analyte Calibrators and Controls Set B

<i>Works with 4 EIAs</i>	#0832 Low Calibrator	#0837 Level I Control	#0833 Cutoff Calibrator	#0838 Level II Control	#0834 Intermediate Calibrator	#0835 High Calibrator
MDMA/Ecstasy (MDMA/Ecstasy EIA)	250 ng/mL	375 ng/mL	500 ng/mL	625 ng/mL	750 ng/mL	1000 ng/mL
Morphine (Opiate 2000 EIA)	1000 ng/mL	1500 ng/mL	2000 ng/mL	2500 ng/mL	4000 ng/mL	6000 ng/mL
Oxazepam (Benzodiazepine 200 EIA)	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
Secobarbital (Barbiturate 200 EIA)	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL

Multi-Analyte Calibrators and Controls Set C

<i>Works with 7 EIAs</i>	#0852 Low Calibrator	#0857 Level I Control	#0853 Cutoff Calibrator	#0858 Level II Control	#0854 Intermediate Calibrator	#0855 High Calibrator
Benzoyllecgonine (Cocaine 150 EIA)	75 ng/mL	112.5 ng/mL	150 ng/mL	187.5 ng/mL	300 ng/mL	1000 ng/mL
d-Methamphetamine (Amphetamines 500 EIA & Methamphetamine EIA)	250 ng/mL	375 ng/mL	500 ng/mL	625 ng/mL	750 ng/mL	1000 ng/mL
Methadone (Methadone EIA)	150 ng/mL	225 ng/mL	300 ng/mL	375 ng/mL	600 ng/mL	1000 ng/mL
Morphine (Opiate 2000 EIA)	1000 ng/mL	1500 ng/mL	2000 ng/mL	2500 ng/mL	4000 ng/mL	6000 ng/mL
Oxazepam (Benzodiazepine 200 EIA)	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL
Secobarbital (Barbiturate 200 EIA)	100 ng/mL	150 ng/mL	200 ng/mL	250 ng/mL	500 ng/mL	1000 ng/mL